

Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023

Page: 1 von 10

SECTION 1: Identification of the substance/preparation and of the company/undertaking**1.1. Product identifier**Name of product: **IMPRIMO® LC MJF**

UFI:

1.2. Relevant identified uses of the substance or mixture and uses advised against**Use of the substance/mixture**

Methacrylate based resin for DLP-systems with 385 nm light sources.

Uses advised against

Private households (= general public)

1.3. Details of the supplier of the safety data sheet**Manufacturer / distributor**

SCHEU-DENTAL GmbH

Am Burgberg 20

D-58642 Iserlohn

Tel. +49 (2374) 9288-0

Fax +49 (2374) 9288-90

eMail: service@SCHEU-DENTAL.comInternet: www.SCHEU-DENTAL.com**1.4. Emergency telephone number:** +49 (2374) 9288-0**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Classification according to Regulation (EC) No. 1272/2008 [CLP]**

Skin Irrit. 2, H315

Eye Irrit. 2, H319

Skin Sens. 1/1A/1B, H317

STOT SE 3, H335

2.2. Label elements

Signal word: Warning

Pictograms: GHS07

**Hazard statements**

H315	Causes skin irritation
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation
H335	May cause respiratory irritation.

Precautionary statements

P261	Avoid breathing dust/fume/gas/mist/vapours/spray
P264	Wash thoroughly after handling.
P271	Use only outdoors or in a well-ventilated area.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P302+352	IF ON SKIN: Wash with plenty of water/soap.
P304+340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+351+338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Material Safety Data Sheet according to (EC) 1907/2006



Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023

Page: 2 von 10

P312	Call a POISON CENTER/doctor if you feel unwell.
P321	Specific treatment (see label).
P332+313	If skin irritation occurs: Get medical advice/attention.
P333+313	If skin irritation or rash occurs: Get medical advice/attention.
P337+313	If eye irritation persists: Get medical advice/attention.
P362+364	Take off contaminated clothing and wash it before reuse.
P403+233	Store in a well-ventilated place. Keep container tightly closed.
P405	Store locked up.
P501	Dispose of contents/container acc. to official regulations.

Special labelling of particular preparations:

none

2.3. Other hazards

No risks worthy of mention.

SECTION 3: Composition / information on ingredients

3.1. Substances

not applicable

3.2. Mixtures

Mixture with, among others, the following ingredients and other non-hazardous admixtures

Composition/information on ingredients

Substance:	CAS-No.:	REACH-No.:	Concentration:	Classification: EC 1272/2008 (CLP)	
Esterification products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enoic acid	41637-38-1		> 95	Skin Irrit. 2, H315; Skin Sens. 1/1A/1B, H317; Eye Irrit. 2, H319; STOT SE 3, H335	ATE (dermal) = Kein Wert ermittelbar; ATE (oral) = Kein Wert ermittelbar; ATE (inhalativ) = Kein Wert ermittelbar
Diphenyl(2,4,6-trimethylbenzoyl)phosphin oxid	75980-60-8		< 2	Repr. 2, H361	ATE (dermal) = Kein Wert ermittelbar; ATE (oral) = Kein Wert ermittelbar; ATE (inhalativ) = Kein Wert ermittelbar

Full text of H- and EUH-phrases: see Section 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

General Information: In case of accident or unwellness, seek medical advice immediately. Take off immediately all contaminated clothing.

In case of inhalation: Provide fresh air. Seek medical attention if problems persist. Seek medical attention if problems persist.

Following skin contact: After contact with skin, wash immediately with plenty of water and soap. In case of skin irritation, consult a physician.

After eyes contact: In case of contact with eyes rinse immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart. Seek medical attention if problems persist.

After ingestion: Rinse mouth immediately and drink plenty of water. No administration in cases of unconsciousness or cramps. Do NOT induce vomiting. Seek medical attention if problems persist.

4.2. Most important symptoms and effects, both acute and delayed

May cause sensitization by skin contact.

Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023
Page: 3 von 10

4.3. Indication of any immediate medical attention and special treatment needed

First Aid, decontamination, treatment of symptoms.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Carbon dioxide (CO₂). Dry extinguishing powder. Water spray jet. Sand

Unsuitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

In case of fire may be liberated: Carbon dioxide (CO₂). Carbon monoxide.

5.3. Advice for firefighters

General information

Co-ordinate fire-fighting measures to the fire surroundings.

Special protective equipment for fire-fighters:

In case of fire: Wear self-contained breathing apparatus. Wear chemical resistant suit.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Wear personal protection equipment. Remove persons to safety. Provide adequate ventilation. Keep away from unprotected people. Keep upwind.

6.2. Environmental precautions

Do not allow to enter into surface or drains. Do not allow to enter into soil/subsoil.

6.3. Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

SECTION 7: Handling and Storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid contact with skin and eyes. Do not breathe gas/fumes/vapour/spray.

Precautions against fire and explosion

Take precautionary measures against static discharges.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed and in a well-ventilated place. Keep container tightly closed in a cool, well-ventilated place. Protect against: UV-radiation/sunlight.

Hints on joint storage

Do not store together with: Oxidizing agent. Keep away from food, drink and animal feedingstuffs.

7.3 Specific end use(s)

Observe instructions for use.

Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023

Page: 4 von 10

SECTION 8: Exposure controls / personal protection
8.1. Control parameters
occupational exposure limit value

Substance	CAS-No.	Source	Occupational exposure limit value: [ppm]	Occupational exposure limit value: [mg/m ³]	Limitation of exposure peaks:	Remark:

Substance with a common (EC) occupational exposure limit value

Substance	CAS-No.	Source	Occupational exposure limit value: [ppm]	Occupational exposure limit value: [mg/m ³]	Limitation of exposure peaks:	Remark:

DNEL/PNEC values
DNEL value

Substance:	CAS-No.:	DNEL/DMEL
Esterification products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enoic acid	41637-38-1	worker inhalative short-term, systemic 3,52 mg/m ³ worker dermal short-term, systemic 2 mg/m ³ /24h population inhalative short-term, systemic 870 µg/m ³ population dermal short-term, systemic 1 mg/m ³ /24h population short-term, systemic 500 µg/m ³ /24h

PNEC-Wert

Substance:	CAS-No.:	PNEC

Additional information

None.

8.2. Exposure controls
Occupational exposure controls.

Provide adequate ventilation as well as local exhaustion at critical locations. Technical measures and the application of suitable work processes have priority over personal protection equipment.

General protection and hygiene measures

When using do not eat, drink, smoke, sniff. Wash hands before breaks and after work. Apply skin care products after work. Wash contaminated clothing prior to re-use.

Personal protection equipment

Minimum standard for preventive measures while handling with working materials are specified in the TRGS 500.

Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required. Respiratory protection necessary at: insufficient ventilation. Suitable respiratory protection apparatus: Protective respiration apparatus not using surrounding air (breathing apparatus) (DIN EN 133).

Hand protection

Tested protective gloves are to be worn: DIN-/EN-Norms: EN ISO 374 Suitable material: Butyl rubber.

Eye/face protection

Eye glasses with side protection

Body protection:

Lab coat

Environmental exposure controls.

See section 7. No additional measures necessary.

Consumer exposure controls

Refer to chapter 7. No further action is necessary.

Exposure Scenario

None.

Material Safety Data Sheet according to (EC) 1907/2006

Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023

Page: 5 von 10

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties****Appearance**

Physical state:	Liquid
Colour:	Transparent
Odour:	Characteristic
Odour threshold:	-

Safety relevant basis data

	Parameter	Value	Unit	Remark
Melting point/freezing point:				No data available
Initial boiling point and boiling range:				No data available
Flammability:				No data available
Lower flammability or explosive limits:				No data available
Upper flammability or explosive limits:				No data available
Decomposition temperature:				No data available
pH:				No data available
Kinematic viscosity:		750	mPa*s	23 °C, 10/s
Water solubility (g/l):				
Partition coefficient: n-octanol/water:				No data available
Vapour pressure:				No data available
Density:		1,1	g/cm ³	
Relative density:				No data available
Particle properties:				No data available

9.2. Other information

None.

SECTION 10: Stability and reactivity**10.1. Reactivity**

No data available.

10.2. Chemical stability

With proper storage and handling the product is stable.

10.3. Possibility of hazardous reactions

No hazardous reaction when handled and stored according to provisions.

10.4. Conditions to avoid

Heat. UV-radiation/sunlight

10.5. Incompatible materials

Can polymerise exothermically if heated, exposed to air, sunlight or by addition of free radical initiators.

10.6. Hazardous decomposition products

In case of fire may be liberated. Pyrolysis products, toxic.

SECTION 11: Toxicological Information**11.1. Information on toxicological effects**

There are no data available on the preparation/mixture itself.

M-factor: -**Acute toxicity (oral):** -**Acute toxicity (dermal):** -**Acute toxicity (inhalative):** -

Material Safety Data Sheet according to (EC) 1907/2006



Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023
Page: 6 von 10

Acute toxicity

Substance:	CAS-No.:	Toxicological information
Esterification products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enoic acid	41637-38-1	LD50 oral (rat) 2000 mg/kg LD50 dermal (rat) 2000 mg/kg

Skin corrosion/irritation:

May cause an allergic skin reaction.

Serious eye damage/irritation:

Strongly irritant. Risk of serious damage to eyes.

Respiratory of skin sensitisation:

May cause sensitization especially in sensitive humans.

CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)

Carcinogenicity:

No indications of human carcinogenicity exist.

Germ cell mutagenicity:

No indications of human germ cell mutagenicity exist.

Reproductive toxicity:

No indications of human reproductive toxicity exist.

STOT-single exposure:

No information available.

STOT-repeated exposure:

No information available.

Aspiration hazard:

No information available.

SECTION 12: Ecological Information

12.1. Toxicity

There are no data available on the preparation/mixture itself.

Ecotoxicity

Substance:	CAS-No.:	Ecotoxicity
Esterification products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enoic acid	41637-38-1	LC50 oral (fish, 96h) 100 mg/l LC50 (crustaceans, 48 h) 6 mg/l EC50 (algae, 72 h) 100 mg/l

12.2. Persistence and degradability

There are no data available on the mixture itself.

12.3. Bioaccumulative potential

There are no data available on the mixture itself.

12.4. Mobility in soil

There are no data available on the mixture itself.

12.5. Results of PBT and vPvB assessment

There are no data available on the mixture itself.

12.6. Endocrine disruptive effect

There are no data available on the mixture itself.

12.7. Other adverse effects

There are no data available on the mixture itself.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Material Safety Data Sheet according to (EC) 1907/2006



Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023

Page: 7 von 10

Appropriate disposal / Product:

Dispose of waste according to applicable legislation. Waste disposal according to directive 2008/98/EC, covering waste and dangerous waste. Must not be disposed together with household garbage.

Appropriate disposal / Package

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

List of proposed waste codes/waste designations in accordance with EWC / AVV

According to EAKV, allocation of waste identity numbers/waste descriptions must be carried out in a specific way for every industry and process.

SECTION 14: Transport Information

Land transport (ADR/RID)

14.1. UN number:	No hazardous material according to transport regulations.
14.2. UN proper shipping name:	
14.3. Transport hazard class(es):	
14.4. Packing group:	
Hazard label:	
Classification code:	
Special Provisions:	
Limited quantity:	
Transport category:	
Hazard No:	
Tunnel restriction code:	
Other applicable information (land transport)	

Inland waterways transport (ADN)

14.1. UN number:	No hazardous material according to transport regulations.
14.2. UN proper shipping name:	
14.3. Transport hazard class(es):	
14.4. Packing group:	
Hazard label:	
Classification code:	
Spezial Provisions:	
Limited quantity:	
Other applicable information (inland waterways transport)	

Marine transport (IMDG)

14.1. UN number:	No hazardous material according to transport regulations.
14.2. UN proper shipping name:	
14.3. Transport hazard class(es):	
14.4. Packaging group:	
Hazard label:	
Special Provisions::	
Limited quantity:	
EmS:	
Other applicable information (marine transport):	

Air transport (ICAO)

14.1. UN number:	No hazardous material according to transport regulations.
14.2. UN proper shipping name:	
14.3. Transport hazard class(es):	
14.4. Packaging group:	
Hazard label:	
Special Provisions:	
Limited quantity Passenger:	
IATA-packing instructions - Passenger:	
IATA-max. quantity - Passenger:	
IATA-packing instructions - Cargo:	
IATA-max. quantity - Cargo:	

Material Safety Data Sheet according to (EC) 1907/2006



Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023
Page: 8 von 10

Other applicable information (air transport)

Passenger-LQ:

14.5. Environmental hazards

ADR/RID / IMDG / ICAO-TI / IATA-DGR: no

Marine pollutant: no

14.6. Special precautions for user

Land transport (ADR/RID)

Transport category: -

Special provisions: -

Tunnel restriction code: -

Limited quantity (LQ): -

Sea transport (IMDG)

EmS-No: -

Special provisions: -

Limited quantity (LQ): -

14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Remark none

SECTION 15: Regulatory Information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU legislation

Information on Regulation (EC) No 166/2006 establishing a European Pollutant Release and Transfer Register:

-

Regulation (EC) No. 1005/2009 on substances that lead to the depletion of the ozone layer:

-

Regulation (EC) No. 648/2004 (Detergents regulation)

-

Regulation (EC) No. 850/2004 [POP-Regulation]:

-

Regulation (EC) No. 649/2012 on the export and import of dangerous chemicals:

-

Use restriction according to REACH annex XVII, no.:

-

National regulatory information

Observe in addition any national regulations!

Restrictions of occupation

none

Other regulations, restrictions and prohibition regulations

-

15.2. Chemical Safety Assessment

For this preparation a chemical safety assessment has been carried out. -

Chemical safety assessments for substances in this mixture were not carried out.

SECTION 16: Other Information

Relevant H- and EUH-phrases (number and full text)

Hazard statements

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

Material Safety Data Sheet according to (EC) 1907/2006



Tradename:

IMPRIMO® LC MJF

Status: 25.01.2023
Page: 9 von 10

H319 Causes serious eye irritation.
H335 May cause respiratory irritation.

Training advice

none

Recommended restrictions of use:

Refer to chapter 1

Further information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

Key literature references and sources for data

Data arise from reference works and literature.

Abbreviations and acronyms

ADR	Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG	International Maritime Code for Dangerous Goods
IATA	International Air Transport Association
GHS	Globally Harmonized System of Classification and Labelling of Chemicals
EINECS	European Inventory of Existing Commercial Chemical Substances
ELINCS	European List of Notified Chemical Substances
CAS	Chemical Abstracts Service
LC50	Lethal concentration, 50 %
LD50	Lethal dose, 50 %